

Caliber – Automated Calibration SW

Changelog

Version 3.03

- Visual Help update
- Fixed CSV and XML data export truncation

Version 2.54 – 3.02

Major updates:

- CSV and XML data export
- ILAC G8 Decision Rules and Guard Band support
- Uncertainty budget overview and result chart for each cal point
- Online licensing

Other updates:

- CAL file editor for calibration data entry
- Custom shortcut key definitions in Modify > Configuration
- Delay countdown shown in procedure window
- Number of discarded samples is now adjustable (1 by default)
- Increased calculation formula length limit
- Support of non-alphanumeric filenames
- Improved analog meter calibration interface
- Calibration protocol updates:
 - o Language profiles for multilingual protocol export
 - o DUT deviation polarity selection for meters and sources
 - o Temp. & Humidity readout from RS232
 - o Serial number readout from Open macro or manual

entry

- o Increased header & footer note length limit
- o „UUT“ changed to „DUT“